

# TOC 1500

Aire cero sin CO2

Número de pieza: 42-1030




## Descripción

Diseñado para proporcionar un suministro de aire confiable para los analizadores de TOC, el TOC1500HP es un sistema robusto utilizado en laboratorios para producir CO2 limpio, seco y libre de CH4



## Generator Details

- **Tipos de gas:** Aire
- **Compresor integrado:** No



## Características clave

- Produce aire limpio, seco, libre de CO2 y CH4
- Ideal para analizadores NPOC TOC
- Funciona 24 horas al día, 7 días a la semana si es necesario. Montaje en pared para laboratorios pequeños
- Pocas piezas móviles significan poco mantenimiento requerido y aseguran una larga vida del generador
- Funcionamiento completamente silencioso
- Se requiere una configuración mínima
- Garantía integral in situ de 12 meses

## Especificaciones técnicas

### TOC 1500 Generador de carbono orgánico completo

Tipo de gas	Aire
Flujo de gas máximo:	1500cc/min
Presión máxima de salida:	80psi/5,5bar
Presión máxima de salida:	80psi/5,5bar
Pureza máxima:	<1 ppm CO2/ 2.0ppm THC
Presión mínima de entrada:	90psi/6,2bar
Presión máxima de entrada:	145psi/10bar
Requisitos de pureza del agua:	N/A
Consumo de agua:	TBC
Voltaje:	100-230 ± 10% VAC
Frecuencia:	50 / 60 Hz
Corriente:	8 Amps
Acreditaciones:	CE
Tamaño (HxAxD) en mm	660 x 250 x 160 mm
Tamaño (HxAxD) en pulgadas	26 x 9,8 x 6,3 inches
Peso del generador	12kg / 26,4lbs

## [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

